

# Product Bulletin

## FoamDDI 1.16 Software Release

Partners & End Users - All Regions

January, 2025

---

**Summary:** This software update includes 2 new features to improve user experience, and a pre-release of 2 new foam applications for FoamDDI, and a number of other fixes.

### Import/Export Calibrations

Users can now export (and import) the module's calibration to USB for the purpose of back-up or transferring to a different LogicBox.

### New Test Methods

2 new test methods are now available ASTM D1881, and Amines. Please contact Visaya for more information and how to configure the instrument to run these new methods.

### Playback Video Markers

Added markers on the playback video slider to indicate where the sample starts foaming and stops foaming.

### General Improvements

**Enhanced Stability:** Addressed multiple recurring software and hardware connectivity issues to ensure smoother operations.

**Firmware Updates:** Updated camera firmware to eliminate frequent disconnection errors and improve synchronization.

**Diagnostics Screen Upgrades:** Improved the reliability of lighting diagnostics, ensuring consistent performance.

# Fixes

**Timer Accuracy and Realism:** Introduced a percentage rather than a timer for better clarity of test progress.

**Enhanced Connectivity:** Stability improvements to prevent frequent disconnections. This includes Cameras specifically and Modules as a whole.

**Diagnostic Lighting Errors/Issues:** Light switches reworked to remove the “Ghost Clicks” as well as the errors that were caused.

**Module Order Broken:** Module settings were fixed so that changing the module order will now persist without having to restart.

**Module Alias:** Module Alias now completely replaces the top header and when defined, no longer displays the default name.

**False Totalizer & Flow Time-Out Errors:** Previously the totalizer timer was inaccurate, a new timer was implemented that is more accurate and eliminates false errors.

**Seq IV Reporting:** The reporting for Seq IV now properly reports Total Volume Increase, and it has been extended to include the Sample Level (per Method D6082).

**Camera Temperature Error Handling and Display:** Cameras used on existing units might not include a heatsink, and will run hot. In some cases they might overheat. This is mostly relevant during repetitive Seq IV tests. The following changes have been made to address the problem:

- All modules now include cameras with appropriate heat sinks (starting early 2024)
- Thermal management (via software) optimizes cameras for low heat.
- Module warning messages will indicate when the cameras are running hot.

Please contact [Service@ayalytica.com](mailto:Service@ayalytica.com) for more information and a link to download this new revision. This software is still compatible with Gen 1 FoamDDI Modules, although some older modules do not support the audible alert.